

Tracheo clean®

DISPOSABLE TRACHEAL PROTECTION








Easier post-operative care

 **Ceredas**
Your confidence inspires us

Ceredas in close collaboration with several ENT and Pneumology departments, has developed a system of protection for the upper respiratory tract. Intended for tracheostomy and tracheotomy patients, this disposable system not only protects the stoma opening but also serves as a tracheal filter.

Tracheoclean® can be used during the post-surgical period as well as at home after discharge.

A non-woven device, much more practical than a regular dressing Tracheoclean® does not require any special care.



Advantages of TRACHEOCLEAN®

- ✔ Allows for easy and fast placement
- ✔ Protects from dust and all contamination
- ✔ Allows for the filtration, heating and moistening of the air breathed
- ✔ Avoids all breathing difficulties with its integrated tracheal strip
- ✔ Provides non-aggressive contact with the sutures
- ✔ Protects health professionals and simplifies their work
- ✔ Preserves the patient's image

Conditioning:

box of 50 or 100 disposable protections

Composition of the disposable tracheal protection:

70% viscose
30% polyester

