

LARYNGECTOMY TUBE  
AND TRACHEOSTOMA SUPPORT  
FOR PULMONARY AND PHONATORY  
REHABILITATION



CyTubeMajor™

LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

CHGE243-GB - © Ceredas 2009

# CyTubeMajor™

LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

A tube for laryngectomy patients, the CyTubeMajor™ is also a support that, immediately after the surgery or even a long time after surgery, is used to receive the different components of the Cyranose® Global System, in particular the Cyranose® HME\*.

Made of medical silicon, it is extremely comfortable to wear and does not damage the tracheal mucosa of the patient when put in place or removed on a daily basis. Fully corresponding with the movements of the neck, it does not provoke discomfort in the everyday life of the laryngectomy patient.

Available in several lengths and diameters to adapt to all sizes of tracheostoma, the CyTubeMajor™ also comes in a fenestrated version for patients with a voice prosthesis, in order to allow for the passage of air in the prosthesis at the time of the obstruction of the tracheostoma during phonation.

## The + of the CyTubeMajor™

► **During hospitalisation :** The CyTubeMajor™ may be used with :

- The Cyranose®, Heat and Moisture Exchanger.
- The Tracheodouche™\*\*, protection during a shower.
- The CyAdaptor™ with the Cyranose® foam filter.

► **Upon discharge from the hospital after cicatrisation and radiotherapy :**

The CyTubeMajor™ may be used with :

- The Cyranose®, Heat and Moisture Exchanger.
- The CyValveFreedom™, a hands-free valve for phonation, by means of the CyAdaptor™.

► **After the end of radiotherapy :**

For laryngectomy patients that have to continue to wear a tube on an occasional or permanent basis, the CyTubeMajor™ may be used with :

- The Cyranose® Heat and Moisture Exchanger.
- The CyValveFreedom™, a hands-free valve for phonation, by means of the CyAdaptor™.
- The Tracheodouche™\*\*, protection during a shower.

## The main **ADVANTAGES** of the CyTubeMajor™

- **Exceptional suppleness** and true comfort and ease due to the use of silicon.
- **Perfectly adapted** to all types of tracheostoma, in particular the difficult morphologies.
- **Absence of deformation** of the support receiving the elements of the **Cyranose® Global System**.
- **Use perfectly compatible** with radiotherapy sessions.
- **Possible pre-fenestrated** for patients with a voice prosthesis.
- **Immediate pulmonary rehabilitation** by connecting the **Cyranose® HME\***, or phonatory rehabilitation with the **CyValveFreedom™** automatic valve.

\*Heat and Moisture Exchanger

\*\*By checking that the permeability is satisfactory between the tracheostoma and the CyTubeMajor™



### *Immediate pulmonary rehabilitation*

*The CyTubeMajor™, immediately after the operation, provides optimum pulmonary rehabilitation with the Cyranose® HME\**



With the CyAdaptor™ and foam filter





With the CyValveFreedom™ automatic valve



# CyTubeMajor™

LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT



| DESCRIPTION   | REFERENCES |   |  |   |   |
|---|------------|---|--|---|---|
|   | L mm       | Size 8<br>Ø I.D. 9,5 mm<br>Ø O.D. 12 mm | Size 9<br>Ø I.D. 10,5 mm<br>Ø O.D. 13,5 mm | Size 10<br>Ø I.D. 12 mm<br>Ø O.D. 15 mm | Size 12<br>Ø I.D. 13,5 mm<br>Ø O.D. 17 mm |
|  <b>CyTubeMajor™</b><br>Laryngectomy tube and tracheostoma support<br>• Standard     | 27         | CBC0827                                 | CBC0927                                    | CBC1027                                 | CBC1227                                   |
|   | 36         | CBC0836                                 | CBC0936                                    | CBC1036                                 | CBC1236                                   |
|   | 55         | CBC0855                                 | CBC0955                                    | CBC1055                                 | CBC1255                                   |
|  <b>CyTubeMajor™</b><br>Laryngectomy tube and tracheostoma support<br>• Fenestrated | 27         | CBCF0827                                | CBCF0927                                   | CBCF1027                                | CBCF1227                                  |
|   | 36         | CBCF0836                                | CBCF0936                                   | CBCF1036                                | CBCF1236                                  |
|   | 55         | CBCF0855                                | CBCF0955                                   | CBCF1055                                | CBCF1255                                  |



# Cyranose®

## GLOBAL SYSTEM

### GLOBAL SOLUTION FOR PULMONARY AND PHONATORY REHABILITATION AFTER A TOTAL LARYNGECTOMY



**Cyranose®** Global System consists of a series of devices made for practitioners and laryngectomy patients for the pulmonary and phonatory rehabilitation of the latter, with or without a voice prosthesis.

\*Heat and Moisture Exchanger